

3M Science.
Applied to Life.™

3M™ Imprint™ 4 Bite VPS Bite Registration Material

Fast and accurate bite registrations

With its great flow properties and no-slump consistency, Imprint 4 Bite material does not lock into undercuts or interdental areas. Material hardness and low flexibility after setting allows for easy cutting and trimming as well as accurate positioning of stone models.

Saves time with a fast intra-oral setting.

Imprint 4 Bite material offers a convenient working time and a short setting time of only 60 seconds, keeping procedures moving quickly while enhancing patient comfort.



Indications

3M™ Imprint™ 4 Bite Registration material is a universal VPS bite registration material. The material can also be scanned for occlusal recording and registration with CAD/CAM systems.



Scan of bite registration made with Imprint 4 Bite material.


	Mixing and application	Setting time from start of mixing*
3M™ Imprint™ 4 Bite VPS Bite Registration Material	0:20	1:00

*The stated times are valid at 23°C/74°F. Times depend on ambient temperature.


Get more confidence in every step.

- 1** Retract
- 2** Syringe and Mix
- 3** Impress

Preliminary Impression



Bite Registration



Ordering information

Item No.	Product Information
71529	3M™ Imprint™ 4 Bite Refill 2 cartridges – 50ml each, 5 Garant™ Mixing Tips green, 5 Garant™ Intra-oral Tips Bite
71460	3M™ Garant™ Intra-oral Tips Bite Includes 50 Intra-oral Tips Bite



3M™ Imprint™ 4 Bite Refill (71529)



3M Deutschland GmbH
Location Seefeld
3M ESPE · ESPE Platz
82229 Seefeld · Germany
info3mespe@mmm.com
www.3MESPE.com

3M, Garant and Imprint are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Please recycle. © 3M 2016. All rights reserved.