



Science.
Applied to Life.™



Stress less over everyday impressions.

Imprint™ 4 VPS Impression Material

Say “YES” to the VPS impression material with the shortest intra-oral setting time. Move through your day with speed, control and precision. Deliver accurate impressions the first time with the fastest VPS material on the market.

It's simply a better way to make better impressions.



Tray Materials

PRODUCT	DISPENSING SYSTEM	VISCOSITY		SETTING VERSION	MAXIMUM WORKING TIME (23°C/73°F)	MAXIMUM INTRA-ORAL SYRINGING TIME (37°C/98°F)	INTRA-ORAL SETTING TIME (37°C/98°F)
		LOW	HIGH				
Imprint™ 4 Penta™ Putty				Regular Set	1:30	–	2:30
Imprint™ 4 Penta™ Heavy Imprint™ 4 Heavy				Regular Set	2:00	–	2:00
Imprint™ 4 Penta™ Super Quick Heavy Imprint™ 4 Super Quick Heavy				Fast Set	1:15	–	1:15

Wash Materials

Imprint™ 4 Light			Regular Set	–	1:00	2:00
Imprint™ 4 Super Quick Light			Fast Set	–	0:35	1:15
Imprint™ 4 Regular			Regular Set	–	1:00	2:00
Imprint™ 4 Super Quick Regular			Fast Set	–	0:35	1:15



www.3M.com/GreatImpressions

3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA
Phone 1-800-634-2249
Web 3M.com/dental

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
Phone 1-888-363-3685

3M, 3M Science. Applied to Life., Imprint and Penta are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2016. All rights reserved. Garant is not a trademark of 3M.