

Fresh material colors,
easier readability.


Imprint™ 4 VPS Impression Material

NEW!

Available in a wide range of viscosities— from Putty to Light Body—Imprint™ 4 VPS Impression Material offers essential benefits for the 1-step and 2-step impression techniques. All tray materials— including Putty—are available for fast and convenient automatic mixing with the Pentamix™ 3 Automatic Mixing Unit and Pentamix™ 2 Automatic Mixing Unit from 3M ESPE.



Impression made with Imprint™ 4 Penta™ Heavy and Imprint™ 4 Regular VPS Impression Material using the 1-step technique.



1-step technique impression made with Imprint™ 4 Penta™ Super Quick Heavy and Imprint™ 4 Super Quick Light VPS Impression Material.

3M ESPE



NEW!

Imprint™ 4

VPS Impression Material

Speed

With innovative advancements in material science, Imprint™ 4 VPS Impression Material cuts intra-oral setting time nearly in half of its predecessor, Imprint™ 3 VPS Impression Material, and other leading competitive VPS dental impression materials. You can noticeably reduce chair time and stress for patients — without reducing the precision of your impression results.

Control

It's easy to get excited about such a drastic reduction in intra-oral setting time, but what good is speed without control? Imprint 4 materials put you in control with adequate working time — a factor for avoiding problems, and for making good impressions with reliable results.

Precision

Imprint 4 material is super-hydrophilic, and displaces moisture right from the start in the unset state, which helps prevent detail-disrupting voids.

Super hydrophilicity results in

- impressions with precise details that are easy to read
- precise-fitting restorations
- fewer adjustments required on final restorations

Fastest intra-oral setting time on the market!

75 sec.

without compromising working time.



Summary of Advantages

- Shortest intra-oral setting time for real time savings
- Sufficient working time for a stress-free procedure
- Super hydrophilicity for accurate detail reproduction
- Clinically proven tear strength
- Fresh material colors for easier readability
- Patient-pleasing mint flavor

Customer Care Center: 1-800-634-2249 www.3MESPE.com/Imprint4

★★★★★ Review these products at www.3MESPE.com/review

3M ESPE

Dental Products

2510 Conway Avenue
St. Paul, MN 55144-1000 USA
1-800-634-2249

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-888-363-3685

Please recycle. Printed in U.S.A.
70-2013-0593-8

3M, ESPE and Imprint are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2013. All rights reserved.