

3M Science.
Applied to Life.™

3M™ Imprint™ 4 Preliminary VPS Preliminary Impression Material

No dust. No mess.

Now you can eliminate the time-consuming and messy work steps needed for alginate impressions. Imprint™ 4 Preliminary VPS material makes your procedure faster and more convenient.

Just automix.

You can easily dispense with a hand dispenser or 3M™ Pentamix™ Automatic Mixing Unit.



Indications

3M™ Imprint™ 4 Preliminary is designed for all typical alginate indications such as:

- Provisional crown and bridge impressions
- Opposing jaw impressions
- Study model impressions
- Preliminary impressions
- Orthodontic impressions
- Impressions for bleaching trays and mouth guards


Product		Maximum Working Time (23 °C)	Constant Intra-Oral Setting Time (37 °C)
Fast	3M™ Imprint™ 4 Preliminary Super Quick	1:00	1:30
	3M™ Imprint™ 4 Preliminary	1:30	2:00
Regular	3M™ Imprint™ 4 Preliminary	1:30	2:00

A back up for peace of mind.

Impressions show virtually no dimensional change — even after three months. This allows for easy remakes of provisionals if necessary.

A smooth matrix makes provisionals shine.

Get an easy shine every time when your 3M™ Protemp™ Plus provisional is made with Imprint 4 Preliminary material.



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web 3M.com/dental

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada

Phone 1-888-363-3685

Ordering information

Item No.	Product Information
3M™ Penta™ Refills	
71521	3M™ Imprint™ 4 Preliminary Penta™ 2 base pastes – 300 ml each, 2 catalysts – 60 ml each
71522	3M™ Imprint™ 4 Preliminary Penta™ Super Quick 2 base pastes – 300 ml each, 2 catalysts – 60 ml each
3M™ Penta™ Value Packs	
71523	3M™ Imprint™ 4 Preliminary Penta™ 6 base pastes – 300 ml each, 6 catalysts – 60 ml each
71524	3M™ Imprint™ 4 Preliminary Penta™ Super Quick 6 base pastes – 300 ml each, 6 catalysts – 60 ml each
Economy Packs	
71539	3M™ Imprint™ 4 Preliminary Regular 10 – 50ml Cartridges
71540	3M™ Imprint™ 4 Preliminary Super Quick 10 – 50ml Cartridges
Bulk Packs	
71548	3M™ Imprint™ 4 Preliminary Regular 30 – 50ml Cartridges
71549	3M™ Imprint™ 4 Preliminary Super Quick 30 – 50ml Cartridges
3M™ ESPE™ Impression Tray Refills (each containing 10 trays)	
71615	Size S, upper jaw
71616	Size M, upper jaw
71617	Size L, upper jaw
71618	Size S, lower jaw
71619	Size M, lower jaw
71620	Size L, lower jaw
3M™ Pentamix™ Lite Automatic Mixing Unit	
77907	3M™ Pentamix™ Lite Automatic Mixing Unit Includes: 1 Penta™ Cartridge
3M™ Pentamix™ 3 Automatic Mixing Unit	
77876	3M™ Pentamix™ 3 Automatic Mixing Unit
Accessories	
76994	Universal 3M™ Penta™ Cartridge for 3M™ Pentamix™ 2 mixing unit
71527	3M™ Imprint™ 4 Preliminary Penta™ Cartridge for 3M™ Pentamix™ 3 mixing unit
71528	3M™ Imprint™ 4 Preliminary Penta™ Super Quick Cartridge for 3M™ Pentamix™ 3 mixing unit
77944	1 3M™ Penta™ Cartridge for 3M™ Pentamix™ Lite mixing unit
77949	50 3M™ Penta™ Mixing Tips – Red

Customer Care Center: 1-800-634-2249
3M.com/greatimpressions

3M, ESPE, Imprint, Penta, Pentamix and Protemp are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Please recycle. Printed in U.S.A. © 3M 2016. All rights reserved. 70-2013-0661-3 (Rev. 8/2016)