3M[™] Medical Specialties

Preliminary Technical Information Bulletin

#1527 Family of Products Comparison Table

| | 1527 | 1527L | 1527LX ³ | 1527ENP | 1527SP |
|--|------------------------------|------------------------------|--|----------------------------|-----------------------------|
| | | | | | |
| Product MSX# | | Discontinued | MSX-6776 ³ | MSX-6774 | MSX-6326 |
| Adhesive Type | Tackified Acrylate* | Tackified Acrylate | Tackified Acrylate | Tackified Acrylate | Tackified Acrylate |
| Backing Style | Square Perforations | Square Perforations | Square ³ Perforations | Square Emboss | Circular Perforations |
| Release Liner Description | No Liner | 83# (6mils/0.15mm) | 83# (6mils/0.15mm) | 63# (4.9mils/ 0.12mm) | 63# (4.9mils/ 0.12mm) |
| Max. Liner Release (gm/25 mm) | N/A | 30 | 30 | 30 | To be determined |
| Air Porosity, Gurley (sec/100cc/sq in) | < 1 – 15 typical | < 1 – 15 typical | Very low | None | < 1 – 15 typical |
| Conformability | High | High | High | Moderate | Moderate |
| Typical Adhesion to Stainless Steel, ² oz/inch (kg/25 mm) | 25 oz/inch (0.7 kg/25 mm) | 32 oz/inch (0.9 kg/25 mm) | 42 oz/inch¹ (1.2 kg/25 mm) | 37 oz/inch (1.1kg/25mm) | oz/inch¹ (kg/25mm) |
| Printable Liner | N/A | Yes | Yes | Yes | Yes |
| Moisture Vapor Transmission Rate (MVTR) | Typical 6200 gm/sqm/day | Typical 6200 gm/sqm/day | Appx. 120 gm/sqm/day ^{1 3} | Appx.50 gm/sqm/day¹ | Typical 6200 gm/sqm/day¹ |
| Roll Width, Maximum | 48 inches (122 cm) | 48 inches (122 cm) | 48 inches (122 cm) | 48 inches (122 cm) | 9.6 inches (24.4 cm) |
| Roll Length, Maximum | 600 yards (549 meters) | 350 yards (320 meters) | 350 yards (320 meters) | 350 yards (320 meters) | 350 yards (320 meters) |
| Shelf Life | 2 years | 2 years | 1 year | 2 years | To be determined |

Current Product Properties – February 2011

¹ Preliminary typical value based on limited number of production lots and subject to change.

² The listed adhesion values differ due to the differences in the manufacturing processes for the various tapes and the point in the process where the sample is collected for testing. The adhesion value may also be influenced by other construction variables, such as the presence or absence of a liner.

³ Only the backing is perforated

* Modified tackified acrylate adhesive - not identical adhesive to others in this tape family.

Product and Safety Information: User is solely responsible for determining the suitability of 3M samples and products for the intended use including any necessary safety or toxicity assessment. 3M will provide Material Safety Data Sheets and results of toxicity testing upon request. In every case before using any product in full scale production users should conduct their own tests to determine to their own satisfaction whether the product is of acceptable quality and is suitable for their particular purposes under their own operating conditions.

Notice: Nothing contained herein shall be construed to imply the nonexistence of any relevant patents or to constitute a permission, inducement or recommendations to practice any invention covered by any patent, without authority from the owners of this patent.

Warranty Information

All statements, technical information and recommendations herein are based on tests 3M believes to be reliable, but the accuracy or completeness thereof is not guaranteed. 3M warrants only that products will meet 3M's specifications at the time of shipment to the customer. 3M does not offer any other warranty and does not warrant the performance, safety or such other characteristics of Products in combination with other materials. 3M specifically DOES NOT warrant Products for any intended or unintended uses (whether or not foreseeable); for compatibility or suitability with other components or compatibility with any methods of manufacture or conversion. The foregoing warranty is made in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability, fitness for a particular purpose and freedom from non-infringement.

Limitation of Remedies: If products are proven not to meet 3M's specifications, the sole and exclusive remedy available and 3M's only obligation shall be, at 3M's option, to replace such quantity of Products which are proven out of specification or to refund the purchase price paid for Products.

Limitations of Liabilities: The remedies provided herein are exclusive remedies against 3M for any alleged or actual nonconformance to specifications or defect or other failure in products or for 3M's performance of its supply obligations. Under no circumstances is 3M liable for any direct, indirect, incidental, special or consequential damages (including lost profits) in any way related to the product under any theory of law including, but not limited to, negligence and strict liability.

Ordering Information

Call us at Customer Service to place an order: 800-742-1994 (U.S.). Visit our website: www.3M.com/medicalspecialties for product and services information, news about conferences we will be attending, new product highlights or to make a direct inquiry.

To have a sales representative contact you, to request samples or clinical and safety summaries, please contact us at 3M HELPLINE 800-228-3957 (U.S.) or for international inquiries, please contact your local country representative. Our 3M Medical Specialties subsidiary contacts are listed below.

International Locations

| Argentina | Dominican Republic | Japan | Poland | Thailand |
|-----------------------------------|---------------------------------|-----------------------------|--|--|
| 3M Argentina S.A.C.I.F.I.A. | 3M Dominicana | 3M Health Care Ltd. | 3M Poland Sp.z o.o. | 3M Thailand Ltd. |
| Tel.: 54-11-4339-2400 | Tel: (809) 530-6560 X499 | Tel: 81-3-3709-8289 | Tel: (48 22) 739 60 00 | Tel: (66)2 260 8577 Ext.175 |
| Fax: 54-11-4339-2640 | Fax: (809) 530-2960 | Fax: 81-3-3709-8754 | Fax: (48 22) 739 60 01 | Fax: (66)2 261 7535 |
| Australia | Finland | Kazakhstan | Portugal | Taiwan |
| 3M Health Care | Suomen 3M Oy | 3M Representation Office | 3M Portugal | 3M Taiwan Ltd. |
| Tel: 1 300 363 878 | Tel.: 358-9-52-521 | Tel: +7 3272 509 944 | Tel: 351-21-3134500 | Tel: 886-2-2704-9011 |
| Fax: 1 800 060 888 | Fax: 358-9-512-29-44 | Fax: +7 3272 509 573 | Fax: 351-21-3134680 | Fax: 886-2-2706-0355 |
| Belgium | France | Korea | Russia | Turkey |
| 3M Belgium N.V./S.A. | Laboratoires 3M Santé | 3M Korea Ltd. | 3M Russia | 3M Turkey |
| Tel: 32-2-722-5111 | Tel: 33-1-3031-8376 | Tel: 82-2-3771-4286 | Tel: 7 (495) 784-74-74 | Tel: 9 (0212) 3507777 |
| Fax: 32-2-720-0225 | Fax: 33-1-3031-8378 | Fax: 82-2-786-2825 | Fax: 7 (495) 784-74-75 | Fax: 9 (0212) 2821742 |
| Brazil | Germany | Malaysia | Singapore | Ukraine |
| 3M do Brasil Ltda. | 3M Medica | 3M Malaysia Sdn. Bhd. | 3M Technologies (S) Pte. Ltd. | 3M Ukraine |
| Tel: 55-19-3838-6661 | Tel: 49 2131 14 40 00 | Tel: 6-03-7806 2888 | Tel: 65-6450-8866 | Tel: +380 44 490 57 77 |
| Fax: 55-19-3838-7449 | Fax: 49 2131 14 49 99 | Fax: 6-03-7806 2902 | Fax: 65-6455-2130 | Fax: +380 44 490 57 75 |
| Canada | Hong Kong | Mexico | Sri Lanka | United Kingdom & Ireland |
| 3M Canada | 3M Hong Kong Ltd. | 3M México, S.A. de C.V. | 3M Lanka (Pvt.) Ltd. | 3M Health Care Ltd. |
| Tel: 800 364-3577 | Tel: (852) 2806 6111 | Tel: +52-55-5270-0400 | Tel: (94) 11 2785701 | Tel: 44-1509-611-611 |
| Fax: 800 341-4630 | Fax: (852) 2234 6044 | Fax: +52-55-5270-0433 | Fax: (94) 11 4410083-4 | Fax: 44-1509-237-288 |
| Chile | Hungary | Middle East | Spain | Venezuela |
| 3M Chile S.A. | 3M Hungária Kft. | 3M Gulf Ltd. | 3M España, S.A. | 3M Manufacturera Venezuela |
| Tel: (56-2) 4103000 | Tel: (36-1) 270-7777 | Tel: 971-4-3670-777 | Tel: 34-91-321-6000 | Tel: 58-241-8391911 |
| Fax: (56-2) 4103400 | Fax: (36-1) 320-0951 | Fax: 971-4-3670-699 | Fax: 34-91-321-6002 | Fax: 58-241-8391928 |
| China | India | The Netherlands | Sweden | Viet Nam |
| 3M China Ltd. | 3M India Ltd. | 3M Nederland B.V. | 3M Svenska AB | 3M Viet Nam Ltd. |
| Tel: 86-21-6275-3535 | Tel: +91-80-2223 1414 | Tel: 31-715-450-450 | Tel: +46-8-92 21 00 | Tel: 84 8 640986 |
| Fax: 86-21-5208-2205 | Fax: +91-80-22231431 | Fax: 31-715-450-212 | Fax: +46-8-92 24 22 | Fax: 84 8 641024 |
| Colombia | Indonesia | New Zealand | Switzerland & Eastern Europ | |
| 3M Colombia S.A. | PT 3M Indonesia | 3M New Zealand | 3M (Schweiz) AG | |
| Tel: 57-1 4101666 | Tel: (62-21) 520 3401 | Tel: 64-9-477 4040 | Tel: +41 44 724 92 44 | |
| Fax: 57-1 4161677 | Fax: (62-21) 520 3106 | Fax: 64-9-477 6692 | Fax: +41 44 724 94 80 | |
| Czech Republic 3M Cesko s.r.o. | Israel 3M Israel Ltd. | Peru 3M Peru S.A. | | ISO 9001 ISO 13485 EN 46001 CERTIFIED |

3M Cesko s.r.o. Tel: +420-261 380 111 Fax: +420-261 380 110

Denmark 3Ma/sTel: 45-43-48-0100 Fax: 45-43-96-8596 3M Israel Ltd. Tel: 972-9-9615000 Fax: 972-9-9615050

Italv 3M Italia SPA Tel: +39 2 7035-2406 Fax: +39 2 7035-2484

3M Peru S.A. Tel: 511-224-2728 Fax: 511-224-3171

Philippines **3M Philippines** Tel: (632) 8133781 to 95 Fax: (632) 8145873 to 74



3M Center, Building 275-5W-05 St. Paul, MN 55144-1000 USA Tel.: 800 228-3957



3M Medical Specialties