


3M™ Pentamix™ Impression Material Delivery Improves Office Infection Control

Once the cartridge and mix tip are installed, the 3M™ Pentamix™ Automatic Mixing Units allow impression material to be hygienically dispensed – with just the push of a button.

Hand Mix Procedure Overview



Pentamix delivery eliminates the need for accessories, like spoons and knives, that require reprocessing steps and can introduce contamination into your procedure.

Dispensing Material with a Pentamix Device



A new 3M™ Penta™ Mixing Tip – Red is applied before every procedure.

More to Love About Automatic Mixing – Efficiency and Accuracy

A recent observation conducted in a group practice showed the following improvements during an experimental period with Pentamix delivery (the practice used 3M™ Pentamix™ Lite Automatic Mixing Unit):

- The incidence rate of the faster one-step technique nearly doubled, saving the practice time*
- Participating doctors experienced a reduced retake rate

Ask your 3M sales representative for more information

* Source: 3M internal data

Ordering Information

3M™ Pentamix™ 3 Automatic Mixing Unit

Item no.	Product information
77876	3M™ Pentamix™ 3 Automatic Mixing Unit Universal Cartridge sold separately
77601	3M™ Pentamix™ 3 Wall Mount Kit
77949	3M™ Penta™ Mixing Tips – Red (50)
Cartridges for 3M™ Pentamix™ 3 Automatic Mixing Unit	
P3792	3M™ Impregum™ Penta™ MB – Purple
P3788	3M™ Impregum™ Penta™ Soft MB – Med Purple
P3794	3M™ Impregum™ Penta™ Soft HB – Dk Purple
P3796	3M™ Impregum™ Penta™ Soft LB – Light Pink
P3776	3M™ Impregum™ Penta™ Soft QS MB – Med Purple
P3774	3M™ Impregum™ Penta™ Soft QS HB – Dk Purple
77784	3M™ Imprint™ 4 Penta™ HB – Sky Blue
77785	3M™ Imprint™ 4 Penta™ QS HB – Mustard Yellow
77786	3M™ Imprint™ 4 Penta™ Putty – Teal Green
71545	Universal Cartridge for 3M™ Pentamix™ 3

3M™ Pentamix™ Lite Automatic Mixing Unit

Item no.	Product information
77907	3M™ Pentamix™ Lite Automatic Mixing Unit includes 1 Penta™ Cartridge for Pentamix Lite
77944	Penta™ Cartridge for 3M™ Pentamix™ Lite
77949	3M™ Penta™ Mixing Tips – Red (50)
71210	3M™ Penta™ Elastomer Syringe Kit includes 1 elastomer syringe, 4 syringe barrels, 5 intra-oral tips

The group practice has continued to use 3M™ Impression Material and Pentamix delivery.

“ Using the Pentamix device allows me to save time and get a more accurate impression...the mixing quality is better, and easier. ”

– Dr. Arie Greenberg
Owner, Greenberg Clinics



Warranty

3M warrants this product will be free from defects in material and manufacture at the time of purchase. 3M MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M's sole obligation shall be repair or replacement of the 3M product.

Limitation of Liability

Except where prohibited by law, 3M will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

3M, Impregum, Imprint, Penta and Pentamix are trademarks of 3M or 3M Deutschland GmbH. Please recycle. Printed in U.S.A. © 3M 2021. All rights reserved.



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web 3M.com/dental