

# Procedure Guide

## Sof-Lex™ Finishing and Polishing System

**A full line of finishing and polishing solutions ... from beginning to end.**

The Sof-Lex system is color-coded from dark (coarse) to light (superfine) for an easily followed step-by-step process. Easy, pop-on mounting lets you quickly change discs or spirals—so you can move efficiently through your finishing and polishing procedure.

Gross Reduction

High-Gloss Polishing



**Sof-Lex™ Diamond Polishing System:** Forget the messy paste. Our pre-polishing spiral prepares the restoration for final gloss, while our diamond-impregnated polishing spiral gives your restorations that gorgeous paste-like gloss.

Pre-Polishing



PP

Diamond Polishing



HGP

Gross Reduction

High-Gloss Polishing



**Sof-Lex™ Discs:** Flexible discs conform to contours of teeth

Gross Reduction



C

Contouring



M

Finishing



F

Polishing



SF

Diamond Polishing



HGP



**Sof-Lex™ XT Discs:** Extra-thin discs for easier access to proximal surfaces

Gross Reduction



C

Contouring



M

Finishing



F

Polishing



SF


Gross Reduction

High-Gloss Polishing

**Sof-Lex™ Strips:**

Different grits on each end with uncoated center for easy interproximal insertion


Gross Reduction & Contouring



C

M

Finishing & Polishing



F

SF

# Sof-Lex™ Finishing and Polishing System







Sof-Lex™ Discs



Sof-Lex™ Diamond  
Polishing System



Sof-Lex™ Strips



[3M.com/SofLex](http://3M.com/SofLex)

**3M**

3M Oral Care  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA  
Phone 1-800-634-2249  
Web [3m.com/SofLex](http://3m.com/SofLex)

3M Canada  
Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
Phone 1-888-363-3685

3M, ESPE and Sof-Lex are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. © 3M 2015. All rights reserved. 70-2013-0721-5