



3M™ Attest™ Steam Biological Indicator (BI) and BI PCD

	Product	Cat. No.	Packaging	Indications for Use	Time to Result	Incubator/Auto-reader
Super Rapid Readout		1491 Blue Cap	50 each/box 4 boxes/case	Gravity: 270°F (132°C) = 3 or 10 min. 275°F (135°C) = 3 or 10 min.		 490/490M Auto-readers
		1492V Brown Cap	50 each/box 4 boxes/case	Dynamic-air-removal (Pre-vacuum and SFPP): 270°F (132°C) = 3 or 4 min. 275°F (135°C) = 3 min.		
		1496V	24 packs/box 24 controls/box	Dynamic-air-removal (Pre-vacuum and SFPP): 270°F (132°C) = 4 min. 275°F (135°C) = 3 min.		
		1496VF	24 packs/box 5 controls/box			
		41482V	24 packs/box 24 controls/box	Dynamic-air-removal (Pre-vacuum and SFPP): 270°F (132°C) = 4 min. 275°F (135°C) = 3 min.		
	41482VF	24 packs/box 5 controls/box				
Rapid Readout		1291	50 each/box 4 boxes/case	Gravity: 270°F (132°C) ≥ 3 min.		 290*/390 Auto-readers
		1292	50 each/box 4 boxes/case	Gravity: 250°F (121°C) ≥ 20 min. Vacuum-assisted: 270°F (132°C) ≥ 4 min.		
		1296	25 packs/box 25 controls/box	Gravity: 250°F (121°C) ≥ 40 min. Vacuum-assisted: 270°F (132°C) ≥ 4 min.		
		1296F	25 packs/box 5 controls/box			
		41382	16 packs/box 16 controls/box	Gravity: 250°F (121°C), 30 min. Dynamic-air-removal 270°F (132°C), 4 min.		
	41382F	16 packs/box 4 controls/box				
Standard		1261	100 each/box 4 boxes/case	Gravity: 270°F (132°C)		 116 Incubator
		1261P	25 each/box 4 boxes/case			
		1262	100 each/box 4 boxes/case	Gravity: 250°F (121°C) Vacuum-assisted: 270°F (132°C) Gravity (table-top sterilizers only): 270°F (132°C) > 3 min.		
		1262P	25 each/box 4 boxes/case			
		1276	25 packs/box 25 controls/box	Gravity: 250°F (121°C) ≥ 30 min. Vacuum-assisted: 270°F (132°C) ≥ 4 min.		
		1276F	25 packs/box 5 controls/box			

*3M™ Attest™ Auto-reader 290 is discontinued.

Refer to the product Instructions for Use (IFU) for a complete list of instructions, warnings, and cautions.

3M Health Care
Medical Solutions Division
2510 Conway Avenue N, Bldg 275-4E-01
St. Paul, MN 55144 USA

© 3M 2020. All rights reserved. 3M and Attest are marks and/or registered marks of 3M. Unauthorized use prohibited. 70-2011-8096-8_US Issued Dec, 2020

Questions? Contact our 3M Sterilization Techline 800.441.1922, Option 2

Need help selecting a BI/BI PCD solution for your steam sterilizer(s)? Try our on-line selection aid: go.3M.com/BIselector

