



Fast, easy composite warming



Solventum™ Filtek™ Composite Warmer

Better flow at the touch of a button

Are you ready to bring composite warming into your practice? With a compact and lightweight design, Filtek Composite Warmer is a fast, easy way to help your procedure flow more smoothly.

- **Fast:** warmed material in as little as 2 minutes after a fast start of only 10 minutes
- **Easy:** one-button operation, easy to clean and disinfect
- **Portable:** removable holder maintains heat and can be moved and used at point of care

Initial warming time:	10 min
Capsule warming time:	2 min
Flowable syringe warming time:	5 min
Footprint:	3 in. (76 mm)
Weight:	190 g
Portable:	Removable holder maintains heat for use at point of care
Capacity:	6 capsules and 1 flowable syringe

Why should you warm your composite?



Better flow



Improved adaptation



Reduced extrusion force

A family of safe-to-warm composites

Filtek™ Dental Restoratives

From easier handling to reduced extrusion force, composite warming can help your procedure flow smoothly.

Built on a history of clinical success, select Filtek™ Dental Restoratives have been tested and indicated for composite warming.^{1, 2}

- 3M™ Filtek™ Supreme Ultra Universal Restorative
- 3M™ Filtek™ Supreme Flowable Restorative
- 3M™ Filtek™ Universal Restorative
- 3M™ Filtek™ Easy Match Universal Restorative
- 3M™ Filtek™ One Bulk Fill Restorative
- 3M™ Filtek™ Bulk Fill Flowable Restorative



¹ Warm capsules to up to 70°C (158°F) for up to 1 hour
² Warm flowable syringes to up to 70°C (158°F) for up to 1 hour, up to 25 times

Ordering Information

Warmer Kit	Item ID
Solventum™ Filtek™ Composite Warmer Kit (includes warmer with heating base, composite holder and sleeve, fixed USB-C cord, power plug, wall outlet adapter, and quick start guide)	6510K

Replacement Parts	Item ID
Solventum™ Filtek™ Composite Warmer Holder Assembly Replacement	6520H
Solventum™ Filtek™ Composite Warmer Silicone Sleeve Replacement	6520S
Solventum™ Filtek™ Composite Warmer Power Plug	6520P

