














3M™ Filtek™ Dental Restoratives

A full line of restorative solutions to fit your unique needs.

3M™ Filtek™ Dental Restoratives have been used in more than 2 billion restorations worldwide. Our family of Filtek dental restoratives lets you choose just the right option to fit any clinical need—from highly esthetic, to highly simple for an easier day.

Properties	Esthetics & Versatility		Esthetics & Simplicity		Efficiency & Simplicity	
	3M™ Filtek™ Supreme Ultra Universal Restorative	3M™ Filtek™ Supreme Flowable Restorative	3M™ Filtek™ Easy Match Universal Restorative	3M™ Filtek™ Universal Restorative	3M™ Filtek™ One Bulk Fill Restorative	3M™ Filtek™ Bulk Fill Flowable Restorative
 Primary benefits	<ul style="list-style-type: none"> Exceptional esthetics with natural-looking fluorescence Outstanding strength for anterior and posterior 	<ul style="list-style-type: none"> Virtually no bubbles or “run-on” Ergonomic design Selection of cannula tip sizes (20G, 27G) 	<ul style="list-style-type: none"> Quick, easy and intuitive shade selection without sacrificing opacity Simplified single-layer restorations, with no blocker needed Naturally-adaptive opacity 	<ul style="list-style-type: none"> Universal opacity for consistent esthetics Pink Opaquer to mask dark areas Low shrinkage stress 	<ul style="list-style-type: none"> Increased opacity for improved esthetics One-step placement up to 5mm 	<ul style="list-style-type: none"> Virtually no bubbles or “run-on” One-step placement up to 4mm Ergonomic design
 Number of shades	36 (Dentin, Body, Enamel, Translucent)	12	3	8 (plus Extra White and Pink Opaquer)	5	4
 TRUE nanotechnology	Yes	Yes	Yes	Yes	Yes	No
 Handling	Firm	Flows on demand	Firm	Creamy	Sculptable	Adapts easily
 Polish retention	Excellent					N/A*
 Wear resistance	Excellent					N/A*
 Safe for composite warming and injection molding^{1,2}	Capsules	Syringes & Capsules	Capsules	Capsules	Capsules	Syringes & Capsules
						

¹ Warm capsules to up to 70°C (158°F) for up to 1 hour. ² Warm flowable syringes to up to 70°C (158°F) for up to 1 hour, up to 25 times. ³ Internal data.
*Filtek Bulk Fill Flowable's intended use requires allowing at least 2 mm on the occlusal cavosurface for a universal or posterior composite.

For a full list of indications, please refer to the product Instructions for Use.

Solventum Dental Solutions
2510 Conway Avenue East-
Maplewood, MN 55144
USA

Phone 1-800-634-2249
Web 3M.com/dental

www.3M.com/Filtek

© Solventum 2024. Solventum, the S logo and Filtek are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. 70-2013-8232-5