

3M Science.
Applied to Life.™



3M™ Coban™ 2 Two-Layer Compression System and 3M™ Coban™ 2 Lite Two-Layer Compression System

3M has advanced the science of compression therapy by designing materials engineered with Intelligent Compression Dynamics. These materials, which comprise the Coban 2 Compression System and Coban 2 Lite Compression System, create a conformable, inelastic sleeve that supports the patient's muscle movements for effective venous return and reduction of edema.*

*Refer to Instructions for Use.



3M™ Coban™ 2 Two-Layer Compression System and 3M™ Coban™ 2 Lite Two-Layer Compression System

Key benefits



Easy to use: Application is fast and easy learn (For Coban 2 Compression System only)



Stays in place: Provides sustained therapeutic compression without slippage



Comfortable to wear: High patient compliance due to product's comfort during wear and low profile



Not made with natural rubber latex



Conforms to the patient's anatomy: Accommodates a variety of shapes and sizes



Provides consistent compression: Compression layer is designed to be applied at full stretch, reducing application variability



Enables a more normal lifestyle: The thin, lightweight, breathable sleeve allows patients to wear their own shoes, so they can return to their regular daily activities

Suggested applications for compression therapy

- ▶ Management of venous leg ulcers, edema, lymphedema and other clinical conditions where compression is appropriate
- ▶ Coban 2 Compression System (Cat #2094N) may be used for patients with an ABPI greater than or equal to 0.8
- ▶ Coban 2 Lite Compression System (Cat #2794N) may be used for patients who have mixed etiology, an ABPI greater than or equal to 0.5, are new to compression (where tolerance is not known), or are frail or less mobile

Ordering information

Application/Use	Name	Cat. No.	Size	Black Nylon Stocking	Rolls/box Boxes/case	Boxes/Case
Below the knee (ABPI ≥ 0.8)	3M™ Coban™ 2 Two-Layer Compression System	2094N	Comfort Layer: 4 in × 2.9 yd (10 cm × 2,7 m) Compression Layer: 4 in × 5.1 yd (10 cm × 4,7 m) stretched	1	2 rolls (1 of each layer) 8	6
Below the knee (ABPI ≥ 0.5)	3M™ Coban™ 2 Lite Two-Layer Compression System	2794N	Comfort Layer: 4 in × 2.9 yd (10 cm × 2,7 m) Compression Layer: 4 in × 5.1 yd (10 cm × 4,7 m) stretched	1	2 rolls (1 of each layer) 8	6
Below the knee (ABPI ≥ 0.8) (large circumference)	3M™ Coban™ 2 Two-Layer Compression System	2094XL	Comfort Layer: 4 in × 3.8 yd (10 cm × 3,5 m) Compression Layer: 4 in × 6.3 yd (10 cm × 5,8 m) stretched	1	2 rolls (1 of each layer) 8	4
Above the knee	3M™ Coban™ 2 Two-Layer Compression System	20096	Comfort Layer: 6 in × 3.8 yd (15 cm × 3,5 m) Compression Layer: 6 in × 4.9 yd (15 cm × 4,5 m) stretched	—	2 rolls (1 of each layer) 8	4



Contact your local 3M account representative, call the 3M Health Care Helpline at **1-800-364-3577**, or visit **3M.ca/coban2layer** for more information.



Available in Canada from:

3M Canada Company
Medical Solutions Division
 300 Tartan Drive
 London, Ontario N5V 4M9
 Canada
 1-800-364-3577
 3M.ca/Medical

3M Company
 2510 Conway Avenue
 St. Paul, MN 55144
 USA
 1-800-228-3957

© 2023 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited. Used under license in Canada. 2302-25662-E CA_EN_70-2013-1524-2