



Tegaderm™

Wound Care Solutions

## 3M™ Tegaderm™ Silicone Foam Dressing

**Your Challenge:** Silicone foam dressings that don't stay in place as long as you need.

**Solution:** A dressing you'll want to stick with.

3M's most advanced silicone foam adhesive dressing takes on your toughest challenges impacting adhesion and wear time:

- ✓ Innovative spoke delivery system allows easier, accurate application
- ✓ Advanced adhesive allows for "lift and check" interval monitoring
- ✓ Conformable thin, low-profile edge minimizes rolling and lifting
- ✓ Unique multi-layer design with high absorbency prevents pooling and reduces risk of maceration



## Easy to apply. Easy to check.

3M knows adhesives like no one else. 3M™ Tegaderm™ Silicone Foam Dressing stays in place with minimal lifting and rolling, and is an excellent choice for your wound management and pressure injury prevention programs.



### Easy, accurate application

A spoke delivery system enables easy application even in challenging locations like the sacrum. Misapplications can result in costly unplanned dressing changes.



### Advanced, 3M silicone adhesive

A unique, proprietary 3M adhesive allows for “lift and check” monitoring of wounds without loss of adhesive properties. You can lift the dressing and reposition it multiple times.



### Thin, low-profile edge










Highly conformable, the thin low-profile border helps minimize the rolling and lifting that can impact adhesion and wear time\*.



### Absorbent multi-layer design

Our innovative layer technology is designed to absorb and wick moisture away from the skin's surface, helping to minimize the backward fluid migration that can cause maceration, disrupt healing, and prompt additional dressing changes.

Tegaderm Silicone Foam dressing adhesive is gentle on even fragile skin, and does not cause damage to the wound, periwound or surrounding skin.

|   | Product               | Cat Number | HCPCS Code | Size                                   | Dressings/Box | Boxes/Case |
|---|-----------------------|------------|------------|--|---------------|------------|
|    | Non-Bordered Dressing | 90631      | A6210      | 4 in. x 4.25 in.<br>10 cm x 10.8 cm    | 10            | 4          |
|    | Non-Bordered Dressing | 90632      | A6210      | 6 in. x 6 in.<br>15 cm x 15 cm         | 10            | 4          |
|   | Bordered Dressing     | 90643      | A6212      | 2 in. x 2 in.<br>5 cm x 5 cm           | 10            | 6          |
|  | Bordered Dressing     | 90640      | A6212      | 3 in. x 3 in.<br>7.5 cm x 7.5 cm       | 10            | 6          |
|  | Bordered Dressing     | 90641      | A6212      | 4 in. x 4 in.<br>10 cm x 10 cm         | 10            | 6          |
|  | Bordered Dressing     | 90642      | A6213      | 6 in. x 6 in.<br>15 cm x 15 cm         | 10            | 4          |
|  | Heel & Contour        | 90646      | A6212      | 6.5 in. x 6.5 in.<br>16.5 cm x 16.5 cm | 5             | 4          |
|  | Small Sacral          | 90647      | A6213      | 6 in. x 6.75 in.<br>15 cm x 17 cm      | 10            | 4          |
|  | Large Sacral          | 90648      | A6213      | 7.25 in. x 8.75 in.<br>18.5 cm x 22 cm | 5             | 4          |

Find out how our silicone foam solution can help address your challenges.

Contact your local 3M Account Representative, call the 3M Health Care Helpline at 1-800-228-3957, or visit [3M.com/SiliconeFoam](https://www.3M.com/SiliconeFoam) for more information.

Disclaimer: HCPCS Codes have been provided to assist you in the preparation of insurance claims. Please note, however, that the reimbursement information provided by 3M Health Care and its representatives is intended to provide general information relevant to coverage and coding for 3M products. Insurers' reimbursement policies can vary and the use of the codes discussed here does not guarantee that an insurer will cover or pay at any particular level. Health care providers should exercise independent clinical judgment in choosing the codes which most accurately describe the products provided.



3M Company  
2510 Conway Avenue  
St. Paul, MN 55144 USA

Phone 1-800-228-3957  
Web 3M.com/medical

\*3M data on file. LAB-SUPPORT-05-310252.

© 3M 2023. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited. 3M marks used under license in Canada. US\_70-2011-6791-6