

3M Glass Ionomer Solutions

Which is best for your practice?

The unique properties of 3M glass ionomers can fit a wide range of patient situations where simplicity, speed, fluoride release and radiopacity are especially important.

3M[™] Ketac[™] Universal Aplicap[™] Glass Ionomer Restorative, is fast, strong and simple.

To try it for yourself ... ask your 3M representative for a sample and demonstration.

3M[™] Ketac[™] Universal Aplicap[™] Glass Ionomer Restorative

- No need for conditioner or coating
- Bulk fill-no need for layering
- Enables fast, one-step, stress-bearing restorative solutions
- Extended indications for use compared to other glass ionomer materials
- Low stickiness allows an easy placement
- Self adhesive
- Available in 6 shades

3M[™] RotoMix[™] Capsule Mixing Unit, 3M[™] Aplicap[™] Activator, 3M[™] Aplicap[™] Applier



3M[™] Ketac[™] Nano Light-Curing Glass Ionomer Restorative

- The esthetic 3M glass ionomer restorative
- Nano technology inspired by 3M[™] Filtek[™] Supreme Ultra Universal Restorative
- Low wear and high initial polish in comparison with other glass ionomer restoratives
- Significantly higher gloss reflectance after polishing
- Faster, easier dispensing with the innovative Quick Mix Capsule or 3M[™] Clicker[™] Dispenser
- Available in 6 shades

3M[™] Elipar[™] DeepCure-S LED Curing Light, 3M[™] Aplicap[™] Applier



Ordering Information

3M[™] Ketac[™] Universal Aplicap[™] Glass Ionomer Restorative

| Item # | Shade/Product Description | | | | | | |
|--------|--------------------------------|--|--|--|--|--|--|
| | Refills | | | | | | |
| | Includes: 50 Aplicap™ capsules | | | | | | |
| 61085 | White, A2, A3 | | | | | | |
| 61086 | White | | | | | | |
| 61087 | A1 | | | | | | |
| 61088 | A2 | | | | | | |
| 61089 | A3 | | | | | | |
| 61095 | A3.5 | | | | | | |
| 61096 | A4 | | | | | | |

| 61097 | Giant Pack A3 |
|-------|---|
| | Includes: 300 Aplicap [™] capsules |
| 73040 | Aplicap™ Activator |
| 73050 | Aplicap™ Applier |
| 76310 | RotoMix™ Capsule |
| | Mixing Unit |

3M[™] Ketac[™] Nano Light-Curing Glass Ionomer Restorative

| ltem # | Shade/Product Description | ltem # | Shade/Product Description | | | |
|--------------------|---|--------|--|--|--|--|
| | Refills Includes: 20 Quick Mix capsules | 3305TK | Trial Kit—Quick Mix Capsule Kit Includes: | | | |
| 3305A1 3305A2 | A1 A2 | | 10–A2 Quick Mix capsules; 10–A3 Quick Mix capsules; | | | |
| 3305A3 3305A3.5 | A3 A3.5 | | Primer | | | |
| 3305A4 | A4 | | | | | |
| 3305B2 | B2 | | | | | |

Saving time and money is important for every practice.

Which 3M glass ionomer restorative is more economical, covering more indications?

| | Indications | | | | | | | | | | |
|--|-------------|--------------------------|------------|---------|----------|----------------|----------|----------|------------------------|-------------|------------|
| | \Im | $\overline{\mathcal{W}}$ | \bigcirc | | Core | Ŵ | Fissure | Delivery | Conditioner/ Primer | Glaze | |
| Product | Class I | Class II | Class III | Class V | Build-Up | Base/Liner | Sealant | Options | Required | Required | Curing |
| Conventional | | | | | | | | | | | |
| 3M[™] Ketac[™] Universal Aplicap[™] Glass Ionomer Restorative | | | | | | | | • | _ | _ | Self Cure |
| 3M[™] Ketac[™] Molar Quick Aplicap[™] Glass Ionomer Restorative | | • | | | | | | • | Ċ | Optional | Self Cure |
| 3M™ Ketac™ Fil Plus Aplicap™ Glass Ionomer Restorative | | _ | | | _ | _ | | • | Ċ | | Self Cure |
| 3M™ Ketac™ Molar Aplicap™ Glass Ionomer Restorative | | | | | | | | • | Ċ | Optional | Self Cure |
| Resin Modified | | | | | | | | | | | |
| 3M™ Ketac™ Nano Light-Curing Glass Ionomer Restorative | | _ | | | | _ | _ | 1- | P | | Light Cure |
| 3M[™] Photac[™] Fil Quick Aplicap[™] Light-Curing Glass Ionomer Restorative | | | | | _ | _ | | • | _ | | Light Cure |
| 3M™ Vitremer™ Core Buildup/ Restorative | • | • | | | | _ | _ | 18 | P | | Tri Cure |
| Indicated | 1 | Quick Mix Ca | apsule 🕴 | Capsule | 3 | BM™ Clicker™ D | ispenser | Handmix | Conditio | oner/Primer | Glaze |

Refer to Instructions for Use for complete product information.



3M Oral Care 2510 Conway Avenue St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249 Web 3M.com/dental For more information on the product you are most interested in, visit

3M.com/KetacUniversal

3M, Aplicap, Clicker, Ketac, Photac, Vitremer and RotoMix are trademarks of 3M or 3M Deutschland GmbH. All other trademarks herein are the property of their respective owners. Please recycle. Printed in USA. © 3M 2023. All rights reserved. 70-2013-0772-8 (REV 6/23)