

# PREVENA™

INCISION MANAGEMENT SYSTEM

## Protect Surgical Incisions, Minimise Postoperative Risk

PREVENA™ Therapy is intended to manage closed surgical incisions and surrounding tissue in patients at risk for developing post-operative complications, including surgical site infection.\*



**PREVENA™ 125**  
Therapy Unit



**PREVENA PLUS™**  
125 Therapy Unit



**CUSTOMIZABLE™**  
Dressing

**PEEL & PLACE™**  
Dressing – 35cm

**PEEL & PLACE™**  
Dressing – 20cm

**PEEL & PLACE™**  
Dressing – 13cm

# Versatility Meets Security

## PREVENA™ Incision Management System: Indications for Use

PREVENA™ Incision Management System is intended to manage closed surgical incisions and surrounding intact skin in patients at risk for developing postoperative complications, such as infection, by maintaining a closed environment via the application of a negative pressure wound therapy system to the incision. The PREVENA™ Incision Dressing skin interface layer with silver reduces microbial colonisation in the fabric.

- Disposable, easy to use negative pressure therapy system
- Visual/audible alerts for patient safety
- Replaceable canister
- Patient can shower with PREVENA™ Dressing in place

### PREVENA™ Incision Dressing skin interface layer:

- Actively wicks fluid away from the incision
- Contains 0.019% ionic silver, which reduces bacterial colonization in the fabric

## The PREVENA™ Incision Management Therapy System clinical outcomes:

More than 70+ published clinical studies in a variety of surgical incisions and patient populations have demonstrated that the PREVENA™ Incision Management System is associated with positive clinical outcomes compared to standard post-operative dressings.



Surgical Site Infection (SSI) & Seroma

## The PREVENA™ Therapy System manages surgical incisions by:



Helping to hold incision edges together



Removing fluids and infectious materials



Acting as a barrier to external contamination



Delivering continuous negative pressure at -125mmHg up to 7 days



Reducing edema

PREVENA™ Incision Management System Bibliography available upon request.

## Ordering information

Item #	Product Name	Qty
PRE1001	PREVENA™ PEEL & PLACE™ System Kit – 20cm	1
PRE1055	PREVENA™ PEEL & PLACE™ Dressing – 20cm	5
PRE1101	PREVENA™ PEEL & PLACE™ System Kit – 13cm	1
PRE1155	PREVENA™ PEEL & PLACE™ Dressing – 13cm	5
PRE3201	PREVENA PLUS™ PEEL & PLACE™ 35cm System Kit	1
PRE3255	PREVENA PEEL & PLACE™ 35cm Dressings	5

Item #	Product Name	Qty
PRE4001UK	PREVENA PLUS™ CUSTOMIZABLE™ System Kit	1
PRE4055	PREVENA PLUS™ CUSTOMIZABLE™ Dressing	5
PRE4000	PREVENA PLUS™ 125 Therapy Unit	1
PRE1095	PREVENA™ 45ml Canister	5
PRE4095	PREVENA PLUS™ 150ml Canister	5
PRE9090	PREVENA™ V.A.C.® Connector	10

## For more information, contact your local representative

**NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. This material is intended for healthcare professionals.**

© 2020 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited. PRA-PM-SG-00032 (07/20)

