


3M™ ESPE™ Intra-oral Syringe - Purple

Single-Use Syringe for Polyether



It's not just a syringe.
It's precision!

Key features

- Excellent ease of use – easy to load, easy to dispense, ergonomic design
- Less material waste compared to mixing tips used with dispenser gun system
- Precise application leads to excellent accuracy of final restorations
- Very hygienic usage and procedure: no messy clean-up, no disinfection
- Choice of wash material type and brand
- Exact control of needed material amount
- Time efficiency in your daily routine

More precision, less waste, easier handling*: get all in one with 3M ESPE Intra-oral Syringe - Purple.

How can Impregum wash application be improved even further? Here's the answer, and it implies a unique auto-mixing solution: the single-use intra-oral syringe by 3M ESPE. This newly developed device can be preloaded directly from the Garant™ cartridge. It is so slim that it enables smooth and easy access to the entire sulcus. And best of all, it allows for exact individual dosage according to requirements. So be clever, forget about waste and look forward to unprecedented advantages.

The intra-oral syringe features all you can wish for to prepare an accurate and flawless Polyether impression*:

- Economical delivery – exact adjustable dose, no waste of material.
- Very hygienic single-use application – no time consuming messy clean-up, no disinfection needed.
- Optimum time efficiency in your daily routine.
- Choice of wash material – no need to buy extra cartridges to fit the syringe.

3M ESPE

Full marks by the users*


In field tests 3M ESPE Intra-oral Syringe - Purple has already proven its advantages. Users particularly appreciate its ease of use and the considerable saving of material due to waste reduction. As compared to a static mixer, the new syringe reduces the residual amount to a fractional part. Small wonder that it has been rated very good to excellent by an overwhelming majority.

89% give top ratings for less amount of waste in the mixing tip




Source: application test with more than 1.100 applications

88% overall satisfaction with intraoral handling



Source: application test with more than 1.100 applications

84% less waste with Impregum™ Garant™ Soft LB




“With these syringes we have an **easier** and **more convenient** application method for impression materials.”

Dr. M. Koroluk, Poland

“**Very nice to use**, much **easier** to dispense compared to bigger gun systems used normally!”

Dr. M. Gregory, UK



Intra-oral Syringe Ordering Information

Item No.	Product Information
71507	3M ESPE Intra-Oral Syringe - Purple Refill Pack 20 Syringes
71508	3M ESPE Intra-Oral Syringe - Purple Value Pack 50 Syringes

3M Australia Pty Limited.
Building A, 1 Rivett Rd
North Ryde NSW 2113.
ABN: 90 000 100 096
Ph 1300 363 454

3M New Zealand Limited.
94 Apollo Drive Rosedale,
Auckland 0632
Ph 0800 80 81 82

*3M ESPE Data on file.

3M, ESPE, Garant, Impregum, Imprint, Penta, Pentamix and Position are trademarks of 3M or 3M ESPE.

Please Recycle. Printed in Australia© 3M 2014.
All rights reserved. PB6636/0214

3M ESPE