













3M™ RelyX™ Universal Resin Cement

Pre-treatment and priming recommendations by material class and indication¹

		Zirconia, metal	Glass-ceramic	Resin based					
		 Crown, bridge, post	 Inlay, onlay	 Adhesive bridge, veneer					
		 Crown, bridge	 Inlay, onlay	 Veneer, tabletop					
		 Crown	 Inlay, onlay	 Veneer, tabletop					
		 Fiber post							
		Sandblast Alumina 30 – 50 µm, 2 bar (30 PSI)	Etch Hydrofluoric acid	Sandblast Alumina 30 – 50 µm, 2 bar (30 PSI)					
3M™ Scotchbond™ Universal Plus Adhesive  <ul style="list-style-type: none"> • Rub in 20 sec • Gently air dry 5 sec 	1	Pre-treat restoration							
	2	Prime restoration							
	3	Etch tooth Phosphoric acid							
	4	Prime tooth							
	5	Apply cement and seat	3M™ RelyX™ Universal Resin Cement 						

¹Main recommendations for typical cases, other options apply. Before using the products described, please refer to the instructions for use.

²Priming not required for 3M™ RelyX™ Fiber Posts

³Optional for cut enamel

Ordering information

3M™ RelyX™ Universal Resin Cement

Item no. Product information

Intro Kit

56968 4x 3M™ RelyX™ Universal Resin Cement syringes shades TR, A1, A3O, WO (3.4g each), 60x 3M™ RelyX™ Universal Micro Mixing Tips, 20 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (5 ml), 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides

Trial Kits

Each trial kit contains: 1x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (1.5 ml), 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides

56969 **Shade TR** **56970** **Shade A1**

Refills

Each refill contains: 1x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips

56971 **Shade TR** **56973** **Shade A3O**

56972 **Shade A1** **56974** **Shade WO**

Value Packs

Each value pack contains: 3x 3M™ RelyX™ Universal Resin Cement syringes (3.4g each), 45x 3M™ RelyX™ Universal Micro Mixing Tips, 15 Elongation Tips

56977 **Shade TR** **56978** **Shade A1**

Mixing Tips

56975 30x 3M™ RelyX™ Universal Micro Mixing Tips

56976 15x 3M™ RelyX™ Universal Micro Mixing Tips with Elongation Tips



3M™ Scotchbond™ Universal Plus Adhesive

Item no. Product information

Vial

41293 **Intro Kit:** 1x 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml), 50 disposable brushes, 1x 3M™ Scotchbond™ Universal Etchant (3ml), 25 dispensing tips

41294 **Refill:** 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml)
1x 5 ml vial **41295** **3x 5 ml vial**

Unit Dose

3M™ Scotchbond™ Universal Plus Adhesive Unit Dose L-Pops, 1x 3M™ Scotchbond™ Universal Etchant (3ml), 25 dispensing tips

41297 **Intro Kit: 50x L-Pop™** **41299** **200x L-Pop™**

41298 **100x L-Pop™** **41304** **400x L-Pop™**

3M™ Scotchbond™ Universal Etchant Etching Gel

41263 2x 3M™ Scotchbond™ Universal Etchant (3ml each), 50 dispensing tips

