



Ambulatory Surgery Center Checklist






We know how important productivity is in your world. We also know it can't come at the expense of quality care. Solventum's proven solutions, supported by decades of experience supporting ambulatory surgery centers, help you minimize the risk of complications and improve efficiency in your facility. In addition, our support programs have proven to be helpful in training staff and streamlining surgery center operations. Throughout the surgical journey, Solventum is here to help.

Product	Product #
Sterilization Monitoring	
As a leader in sterilization monitoring assurance for decades, it is your needs and challenges that drive us to keep innovating. Our sterilization monitoring solutions are designed to simplify, standardize, and streamline workflows.	
3M™ Attest™ Super Rapid Readout Biological Indicators Super rapid biological indicators give both steam and VH2O2 results in 24 minutes*	1491 (brown cap) 1492V (blue cap) 1295 (pink cap) 24 minute results* 
3M™ Attest™ Auto-readers Intelligent display with time-to-result countdown feature. Display can be viewed remotely using your PC using Attest Connect software..	490M 490
3M™ Attest™ Chemical Indicators Immediate, accurate and easy-to-read monitoring of sterilization process conditions inside packs for use in steam and low temp sterilizers. Distinct indicating readouts for instant results with easy interpretation.	1243A 1243B 1248 1348 
3M™ Steam Indicator Tapes Trusted, high quality adhesive with stretchable backing minimizes tape “pop-off” during sterilization.	1355-18MM 1355-24MM
3M™ Attest™ eBowie-Dick and Bowie Dick Test System Bowie Dick test packs are designed to monitor the performance of dynamic-air-removal steam sterilizers each day. The eBowie-Dick Test System provides automation and digitization, removing guess-work.	10135 (eBowie Dick Test Card) 10135CH (eBowie Dick Card Holder) 1190 (eBowie Dick Auto-Reader) 00135LF (Bowie Dick Plus Test Pack)

Note: Find products and sizes within Sterilization Monitoring brochure or at go.3m.com/attest.

Patient Warming	
Prewarming with Bair Hugger warming blankets and warming gowns is an important step in a comprehensive patient warming protocol that also includes intraoperative warming.	
3M™ Bair Hugger™ Warming Blankets The Bair Hugger warming blanket system offers flexible temperature management solutions designed to help clinicians deliver optimal care while maintaining normothermia throughout the perioperative process.	30000 52500 62200
3M™ Bair Hugger™ Warming Gowns The Bair Hugger warming gown system adds a new dimension to patient warming or prewarming, providing clinical and comfort warming options in a single patient gown that can be used throughout the perioperative process.	81001 81002 81003
3M™ Bair Hugger™ Warming Unit The Bair Hugger warming units are used in conjunction with warming blankets and gowns to help prevent unintended hypothermia.	67500 77500 87500+
3M™ Bair Hugger™ Temperature Monitoring System Non-invasive, accurate, consistent, and continuous core-temperature monitor that can be used throughout the perioperative process and beyond.	36000 37000
3M™ Ranger™ Blood/Fluid Warming System Covers the entire spectrum of fluid warming — from KVO to 30,000 mL/hr — with one warming unit.	24500

* Must be used in conjunction with 3M™ Attest™ Auto-reader 490 or with a Model 490H having software version 4.0.0 or greater or with 3M™ Attest™ Mini Auto-Reader 490M. 24-minute results available only on 490/490H Auto-readers having software versions 4.0.0 or greater.
* 3M™ Attest™ Biological Indicators 1491, 1492V and 1295 can be read in 24 minutes when used with 3M™ Attest™ Mini Auto-reader 490M or with 3M™ Attest™ Auto-reader 490 or 490H with software version 4.0.0 or greater.
+ 87500 is for use in pre-op and post-op only, not for use in OR/intra-op.

Product		Product #
Skin and Nasal Antiseptic		
3M™ Skin and Nasal Antiseptic Take control of preoperative nasal decolonization with 3M™ Skin and Nasal Antiseptic. <ul style="list-style-type: none"> • A pH-balanced formulation with a scientifically developed film-forming polymer to increase persistence • Safe and non-irritating when used in the nares 		192401 
Surgical Clippers		
3M™ Surgical Clipper Professional and 3M™ Surgical Clipper Professional Blade 3M™ Surgical Clipper Professional combines superior design and performance to help reduce the risk of SSIs, streamline workflow and promote a positive perioperative experience for your staff and patients.		9681 Clipper 9680 Blade 
Medical Securement		
<i>Meet securement needs while choosing the right tape for the right patient and application, Solventum has medical tapes for fragile/at-risk skin and intact skin.</i>		
3M™ Medipore™ H Soft Cloth Surgical Tape Soft cloth, conformable tape to maintain skin integrity. Features include multi-direction stretch, perforated and easy to tear.		2860S-2U
3M™ Micropore™ S Surgical Tape* Reliable adhesion and removes with minimal disruption of skin layers and with minimal discomfort.		2770S-1
3M™ Durapore™ Surgical Tape Securement for central venous catheters, foley catheters, nasogastric tubes, chest tubes, orogastric tubes and surgical drain tubes.		1538S-1
3M™ Coban™ NL Non-Latex Self-Adherent Wrap* Used in clinical applications such as blood draws, dressings, immobilization and securement for difficult to dress areas.		2082-1X
		
IV Site Management Solutions		
<i>IV catheter securement and antimicrobial port protection.</i>		
3M™ Curores™ Disinfecting Port Protectors <ul style="list-style-type: none"> • Luer-lock design to provide a brightly colored visible physical barrier to contamination between assessments • Save time, improve consistency and enhance compliance by using passive disinfection of IV access points. 		CM50200 (male luers) CSA5-250 (open female luers)
3M™ Tegaderm™ I.V. Advanced Securement Dressing A breathable, transparent film dressing designed to deliver exceptional patient care through advanced primary securement of IV catheters with enhanced wear time and gentle removal.		1685 (Central lines) 1882 (Peripheral lines) 9132 (with gentle removal and antimicrobial protection)
3M™ Tegaderm™ Transparent Film Dressing with Adhesive-Free Window Tegaderm™ transparent dressings can be used to cover and protect catheter sites and wounds.		1624W 1626W
		
Surgical Hand Antiseptic		
3M™ Avagard™ Surgical and Healthcare Personnel Hand Antiseptic with Moisturizers (Chlorhexidine Gluconate 1% solution and Ethyl Alcohol 61% w/w) A waterless, brushless surgical hand antiseptic. Formulated with two active ingredients — alcohol for rapid kill and chlorhexidine gluconate for persistent activity.		9200 (500ml) 9216 (1200ml) 

*Skin performance: The icon designates products that deliver the securement power you need while minimizing damage to skin.

Product	Product
Surgical Skin Prep Solutions	
<i>There is no one surgical skin prep solution suitable for all surgeries- whether containing iodine, alcohol, CHG or a combination of these. Patient considerations and prep contraindications must be considered. Be prepared with a family of surgical skin preps to help reduce the risk of surgical site infections.</i>	
3M™ SoluPrep™ S Sterile Antiseptic Solution chlorhexidine gluconate [2% w/v] and Isopropyl alcohol [70% v/v] Patient Preoperative Skin Preparation.	67130 (26ml) 67133 (10.5ml)
3M™ DuraPrep™ Surgical Solution (Iodine Povacrylex [0.7% available iodine] and Isopropyl Alcohol, 74% w/w) Patient Preoperative Skin Preparation.	8630 (26ml) 8635 (6ml)



Surgical Drapes	
3M™ Ioban™ 2 Antimicrobial Incise Drapes Reduce bacterial surgical site contamination compared to no incise drape. Provide continuous broad-spectrum antimicrobial activity to help reduce the risk of surgical site contamination that may be associated with surgical site infections.	6640EZ 6648EZ 6650EZ 6651EZ 6661EZ
3M™ Ioban™ CHG Chlorhexidine Gluconate Incise Drape (2% w/w CHG) Building on the heritage of 3M™ Ioban™ 2 Antimicrobial Incise Drapes, this advanced surgical incise drape with chlorhexidine gluconate (CHG) formulated into the adhesive creates a sterile operative surface to help reduce the risk of skin flora contamination.	8840EZ 8850EZ 8851EZ



Electrodes	
3M™ Red Dot™ Electrodes Provide the right electrocardiogram (ECG) electrodes for a variety of patients and monitoring situations.	2670 2570



Incision Management and Dressings	
3M™ Prevena™ Incision Management System Reduces post-operative edema, delivers -125mmHg for therapy length of 7 days with the 3M™ Prevena™ 125 Therapy Unit, helps hold the incision edges together, removes fluids and infectious materials and protects the incision from external infectious sources.	PRE1101US (13cm) PRE1001US (20cm)
3M™ Tegaderm™ Absorbent Clear Acrylic Dressing Offers a waterproof, transparent barrier — allowing for active monitoring without disturbing the incision.	90815
3M™ Tegaderm™ +Pad Dressings with Non-Adherent Pad All-in-one dressing that provides a sterile, waterproof, viral and bacterial barrier. These dressings consist of a non-adherent absorbent pad bonded to a thin film transparent dressing	3582 3586 3591
3M™ Steri-Strip™ Adhesive Skin Closures Adhesive skin closure devices may be used in skin approximation of lacerations and surgical incisions.	R1546 R1547 R1548
3M™ Medipore™ +Pad Soft Cloth Adhesive Dressing Soft, porous and breathable. Unique, non-adherent pad won't stick to the wound, thus minimizing disruption of the healing tissues upon removal. Absorbent pad protects wound and absorbs drainage, while allowing skin to breathe. Helps reduce the risk of skin maceration.	3566 3568 3569 3570 3571



Splinting products	
3M™ Scotchcast™ Quick Step Double Sided Felt Roll Splint Comfortable, customizable and comfortable. Provides strength and stability with pre-padded, breathable felt to ensure quick drying.	74003Q, 3" x 15' (other sizes available)



Hand hygiene	
<i>The right product for every need. Gentle, effective hand hygiene and skin care solutions. 3M™ Avagard™ Hand Antiseptics provide fast-acting solutions to help protect your patients and staff.</i>	
3M™ Avagard™ D Instant Hand Antiseptic with Moisturizers (Ethyl Alcohol 61% w/w)	9221 (88 mL) 9222 (500 mL) 9230 (1000 mL)



Additional Solventum Medical product options and sizes are available. Visit [Solventum.com](https://www.solventum.com) for more information.

Product	Product #
Skin Care	
3M™ Cavilon™ No Sting Barrier Film* The original and only terpolymer-based alcohol-free barrier film** that helps prevent skin damage before it occurs.	3344 3343
3M™ Cavilon™ Moisturizing Hand Lotion A CHG compatible lotion designed to restore moisture and prevent dry skin. Hypoallergenic. Fragrance-free. Compatible with the 3M™ Avagard™ brand of products.	9205 (16 fl oz)



Solventum Education provides free contact hours through online and in-person educational content dedicated to helping clinicians focus on deepening expertise and improving patient care.

For more information visit:
<https://go.solventum.com/education>



Learn more about how Solventum can help you minimize the risk of complications and drive efficiency:

- Contact your local **Solventum Account Manager**
- Call the Solventum Health Care helpline at **800-228-3957**
- Visit **Solventum.com**

**Of leading competitors in the market, based on disclosed ingredient information.
 * Flammable: Consult Instructions for Use for application instructions and dry time criteria.



Solventum Corporation
 2510 Conway Avenue
 St. Paul, MN 55144-1000
 USA

Phone 800-228-3957
Web Solventum.com

Note: Specific instructions, limitations, contraindications, warnings, precautions, and safety information exist for these products and therapies, some of which may be Rx only. Please consult a clinician and product instructions for use prior to application.

© Solventum 2025. Solventum, the S logo, V.A.C., Veraflo and Veraflo Cleanse Choice are trademarks of Solventum or its affiliates. 3M is a trademark of 3M. Other trademarks are the property of their respective owners.

US_70-2011-7154-6